

PRODUCT CODE

BS004

APPLICATION

A Romanowsky stain for blood and bone marrow morphology

REAGENT COMPOSITION

Giemsa Stain Powder	0.7 g
Glycerol	250 ml
Methanol	750 ml

REAGENT PREPARATION

It is recommended to filter all stain before use.

REAGENT STORAGE AND STABILITY

Giemsa stain solution is stable up to the stated expiry date when stored at 15-25° C.
Keep tightly closed to prevent air oxidation

PROCEDURE

1. Fix air dried smears 15 minutes in methanol in a staining jar.
2. Transfer smears directly to a staining jar containing Giemsa Stain freshly diluted 1:10 with pH 6.8 buffer or distilled water. Stain for 10- 15 minutes.
3. Wash in several changes of buffered or distilled water, allow drying.

Staining with May Grunwald's Stain:

1. Fix air dried smears 15 minutes in methanol in a staining jar.
2. Transfer smears directly to a staining jar containing May-Grunwald's stain freshly diluted 1:1 with pH 6.8 buffer or distilled water. Stain for 15 minutes.
3. Transfer smears directly to a staining jar containing Giemsa Stain freshly diluted 1:10 with pH 6.8 buffer or distilled water. Stain for 10- 15 minutes.
4. Wash in several changes of buffered or distilled water, allow drying.















NOTE

This method takes much long than Romanowsky staining techniques but the brilliance of the Stained cells makes the extra time worthwhile particularly for bone marrow smears.

RESULTS

Red cells	-	Pink
Cell nuclei	-	Blue- Purple
Acidophil granules	-	Pink- Red
Basophil granules	-	Blue

SYMBOLS ON LABELS

Symbols	Signify	Symbols	Signify
	Catalogue Number		Pack Size
	Expiry Date		Volume
	Storage Condition		Lot Number
	Instruction for Use		In Vitro Diagnostics
	Manufacturing Date		Manufacturer
	Number of Tests		For Single Use Only
	EC Representative		European conformity

BIBLIOGRAPHY

1. Giemsa, G. (1922): Das Wesen der Giemsa-Farbung, Zentralbl f Bakt; p89, pp99-106.
2. May, R. et Grünwald L. (1909): Über die Farbung von Feutchpreparaten mit meiner Azur-Eosine methode, Deutsche med Xschr, p35, pp1751-1752.

