

PRODUCT CODE

CE005

INTENDED USE

This reagent kit is intended for the “In Vitro “quantitative determination of Potassium in Serum

CLINICAL SIGNIFICANCE

Potassium is an electrolyte that is vital to cell metabolism. It helps transport nutrients into cells and removes waste products out of cells. It is also important in muscle function, helping to transmit messages between nerves and muscles. Elevated potassium levels (hyperkalemia) are often associated with renal failure, dehydration shock or adrenal insufficiency. Decreased potassium levels (hypokalemia) are associated with malnutrition, negative nitrogen balance, gastrointestinal fluid losses and hyperactivity of the adrenal cortex, this test measures the amount of potassium in the blood

PRINCIPLE

Potassium reacts with sodium tetra phenol boron in a specially prepared buffer to form a colloidal suspension. The amount of the turbidity produced is directly proportional to the concentration of potassium in the sample.

REAGENT COMPOSITION

Reagent 1: Potassium Reagent

Reagent 2: Potassium Standard 5 mEq/L

REAGENT PREPARATION

The reagent and standard are ready to use.

STORAGE AND STABILITY

The reagents and standard are stable up to the stated expiry date when stored at 15-30 ° C.

SPECIMEN

Serum free from hemolysis

Separate serum from the clot as soon as possible as Potassium may leach from the RBC's which can elevate potassium level.

NORMAL RANGE

Serum/Plasma: 3.5 - 5.5 mEq/L.

It is recommended that each laboratory establish its own normal range representing its patient population.

ASSAY

Wavelength 630 nm
 reaction type Endpoint
 Cuvette 1 cm light path
 Temperature Room temp.
 Measurement Against reagent blank

PROCEDURE

Pipette in to cuvettes	Blank	Standard	Sample
Potassium Reagent	1000 µL	1000 µL	1000 µL
Standard	--	20 µL	--
Sample	--	--	20 µL

Mix and incubate for 5 min at room temp. And measure the absorbance of the sample (As) and the standard (A.std) against the reagent blank.

CALCULATION

$$\text{Concentration of Potassium (mEq/L)} = \frac{\Delta A \text{ sample}}{\Delta A \text{ standard}} \times 5 (\text{Std. conc.})$$

LINEARITY

This procedure is linear up to 7 mEq/L. If values exceed his limit dilute the sample with distilled water and multiply results with proper dilution factor.















NOTE

As potassium is a very widely distributed ion, care should be taken to avoid any contamination. All glass wares being used for the test should first be rinsed with 1% or 0.1 N HNO3 and then with good quality deionized water before use.

QUALITY CONTROL

Control serum of known concentrations should be analyzed with each run.

SYMBOL ON LABELS

Symbols	Signify	Symbols	Signify
	Catalogue Number		Pack Size
	Expiry Date		Volume
	Storage Condition		Lot Number
	Instruction for Use		In Vitro Diagnostics
	Manufacturing Date		Manufacturer
	Number of Tests		For Single Use Only
	EC Representative		European conformity

BIBLIOGRAPHY

- 1-Tietz, N.W., Fundamentals of Clinical Chemistry, W.b Saunders Co. Phila, P.A. p. 874.
- 2-enry R.F., et, al, Clinical Chemistry Principles and Te2-chnics. 2ndE d, H arper an d R ow, Harper and Row, Hargersein, M.D. (1974)
- 3- Maruna RFL., Clin Chem. Acta. 2:581, (1958)
- 4- Trinder, P:Analyst, 76:596, (1951)