

### PRODUCT CODE

CS001

### INTENDED USE

This reagent is intended for in vitro quantitative determination of Albumin in human serum.

### CLINICAL SIGNIFICANCE

Albumin is the most abundant protein constituent of serum. It is synthesized in the liver and is noted for its ability of configuration changes. This steric affinity allows the albumin molecule to serve as a carrier of many substances such as bilirubin, fatty acids, uric acid, various drugs, and antibiotics. Albumin also functions in the maintenance of proper osmotic pressure.

Elevated serum albumin levels are associated with possible dehydration. Low serum albumin levels are indicative of potential malnutrition, liver disease, kidney disorders, and rheumatoid arthritis.

### PRINCIPLE

Serum albumin binds selectively to the dye bromocresol green at pH 4.2. The absorbance of the blue/green complex at 546 nm is proportional to the albumin concentration.

### REAGENT COMPOSITION

#### Albumin (Liquid) Reagent

Citrate Buffer (pH4.2) 7.5 mmol/L  
Bromocresol Green 150 µmol/L

#### Albumin Standard

Albumin standard concentration 4 g/dL

### REAGENT PREPARATION

Reagent and standard are ready for use.

### REAGENT STORAGE AND STABILITY

The color reagent and standard are stable up to the stated expiry date when stored at 2-8° C. Contamination after opening must be avoided. The reagent should be a clear yellow / green solution. If turbidity or precipitation has occurred the reagent should be discarded.

### SPECIMEN

Serum, heparinized or EDTA plasma  
Stability in serum at 2-8°C for 1 month at 15-25° for 1 week

### PRECAUTION

To avoid contamination, use clean laboratory wares.  
Avoid direct exposure of reagent to light.

### ASSAY

Wavelength : 546nm  
Cuvette : 1 cm light path  
Temperature : 20-25°C  
Measurement : Against reagent blank

### PROCEDURE

Pipette in to cuvette	Blank	Standard	Sample
Albumin Reagent	1000 µL	1000 µL	1000 µL
Standard	--	10 µL	--
Sample	--	--	10 µL

Mix and incubate for 5 minutes at 20-25°C. Measure the absorbance of the sample (As) and standard (Astd) against the reagent blank within 30 minutes.

### CALCULATION

$$\text{Serum Albumin (g/dL)} = \frac{\Delta A \text{ sample}}{\Delta A \text{ standard}} \times 4 \text{ (Std.conc.)}$$

### LINEARITY

The test is linear up to a concentration of 7g/dL. If a higher albumin concentration is expected, dilute sample 1+1 with physiological saline. Repeat the estimation and multiply the result by 2.

### NORMAL RANGE

It is recommended that each laboratory establish its own reference values. The following value may be used as guide line.

Serum Albumin: 3.8 - 5.1 g/dL

### QUALITY CONTROL















All control sera with Albumin value estimated by this method can be used.

### NOTES

- The test is not influence by bilirubin values up to 20 mg/dl.
- In case of excessive lipemic a sample blank should be prepared by adding 0.025 ml of serum to 2.5 ml of 0.9 % saline. The absorbance of the sample blank is subtracted from the absorbance of the sample.
- Avoid excessive haemolysis since every 100 mg/dl of haemoglobin corresponds to about 100 mg/dL of Albumin.

- The color reagent and standard contain sodium azide. Do not swallow. Avoid contact with skin and mucous membrane.

### SYMBOL ON LABELS

Symbols	Signify	Symbols	Signify
	Catalogue Number		Pack Size
	Expiry Date		Volume
	Storage Condition		Lot Number
	Instruction for Use		In Vitro Diagnostics
	Manufacturing Date		Manufacturer
	Number of Tests		For Single Use Only
	EC Representative		European conformity

### BIBLIOGRAPHY

- Doumas . B. T. et al ; Clin. Chim. Acta.. 31.87, 1971.
- Tietz. N.W. (Ed); Text book of clinical Chemistry , W.B. Saunders, 589, 1986
- Doumas . B. T. et al ; Standard methods of Clinical Chemistry , 7, 175 ; Academic Press of Chocago 1972.
- Walsh , R.L.; Clin. Biochem., 16, 178, 1983.

