

**PRODUCT CODE**  
CS025

**INTENDED USE**

The reagent is intended for in vitro quantitative determination of Direct Bilirubin in serum or plasma.

**CLINICAL SIGNIFICANCE**

Bilirubin is a breakdown product of hemoglobin. It is transported from the spleen to the liver and excreted into bile. Hyperbilirubinemia results from the increase of bilirubin concentrations in plasma. Causes of hyperbilirubinemia:  
 Total bilirubin: Increase hemolysis, genetic errors, neonatal jaundice, ineffective erythropoiesis, and drugs.  
 Direct bilirubin: Hepatic cholestasis, genetic errors, hepatocellular damage<sup>1,6,7</sup>  
 Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**PRINCIPLE**

Bilirubin is converted to colored azobilirubin by diazotized sulfanilic acid and measured photometrically. Of the two fractions presents in serum, bilirubin-glucuronide and free bilirubin loosely bound to albumin, only the former reacts directly in aqueous solution (bilirubin direct), while free bilirubin requires solubilization with dimethylsulphoxide (DMSO) to react (bilirubin indirect). In the determination of indirect bilirubin, the direct is also determined; the results correspond to total bilirubin. The intensity of the color formed is proportional to the bilirubin concentration in the sample<sup>1,2,3</sup>

**REAGENT COMPOSITION**

<b>Direct Bilirubin Reagent (R1)</b>	
Sulphanilic Acid	30 mmol/L
Hydrochloric acid (HCl)	150 mmol/L
<b>Direct Bilirubin, Nitrite Reagent (R2)</b>	
Sodium nitrite	29 mmol/L

**PREPARATION**

Both reagents are ready to use.

**REAGENT STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Color development in R 2.

**SPECIMEN**

-Serum or plasma, free of hemolysis  
 -Protect samples from direct light. Stability: Bilirubin is stable at 2-8°C for 4 days and 2 months at -20°C.

**PRECAUTION**

R1: H290-May be corrosive to metals. H314-Causes severe burns and eye damage. EUH208-Contains sulphanilic acid. May produce an allergic reaction.

**PROCEDURE**

**1.ASSAY Condition**

Wavelength	555 nm (530-580)
Temperature	15-25 °C
Measurement	Against sample blank (without nitrite)

**2. Adjust the instrument to zero with distilled water.**

**3.Pipette into a cuvette**

	Blank	Sample
direct Bilirubin Reagent (R1)	1000 µL	1000 µL

direct Bilirubin, Nitrite reagent (R2)	--	50 µL
Sample	100 µL	100 µL

- Mix and stand for exactly '10' minutes at room temperature.
- Measure the absorbance of sample Blank (As).

**CALCULATION**

**-With Calibrator:**

$$\frac{(A)Sample - (A)Sample\ Blnk}{(A)Calibrator - (A)Calibrator\ Blank} \times Conc.\ Calibrator = mg/dl\ of\ bilirubin\ in\ the\ sample$$

**-With Factor :**

$$((A)Sample - (A)Sample\ Blank) \times Factor = mg/dl\ of\ bilirubin\ in\ the\ sample$$

$$Factor = \frac{Concentration\ of\ Calibration}{(A)Calibrator - (A)Calibrator\ Blank}$$

Conversion Fctor : mg/dl x 17.1 µmol.

**NORMAL RANGE**

Bilirubin Direct: Up to 0,25 mg/dL  $\cong$  4,27 µmol/L  
 These values are for orientation purpose; each laboratory should establish its own reference range















**QUALITY CONTROL**

All control sera with Direct Bilirubin value estimated by this method can be used.  
 Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**LINEARITY**

From detection limit of 0,07 mg/dL to linearity limit of 20 mg/dL.  
 If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

**SYMBOL ON LABELS**

Symbols	Signify	Symbols	Signify
	Catalogue Number		Pack Size
	Expiry Date		Volume
	Storage Condition		Lot Number
	Instruction for Use		In Vitro Diagnostics
	Manufacturing Date		Manufacturer
	Number of Tests		For Single Use Only
	EC Representative		European conformity

**BIBLIOGRAPHY**

- David G Levitt and Michael D Levitt. Quantitative assessment of the multiple processes responsible for bilirubin homeostasis in health and disease. Clin Exp Gastroenterol. 2014; 7: 307-328.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.