

PRODUCT CODE

CZ012

INTENDED USE

The reagent is intended for in vitro quantitative determination of creatine kinase in serum or plasma.

CLINICAL SIGNIFICANCE

CK-MB is an enzyme formed by the association of two subunits from muscle (M) and nerve cells (B). CK-MB is usually present in serum at low concentration; it is increased after an acute infarct of myocardium and later descends at normal levels. Also is increased, rarely, in skeletal muscle damage^{5,6,7,8}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

PRINCIPLE

The procedure involves measurement of CK activity in the presence of an antibody to CK-M monomer. This antibody completely inhibits the activity of CKMM and half of the activity of CK-MB while not affecting the B subunit activity of CK-MB and CK-BB. Then it's used the CK method to quantitatively determine CKB activity^{1,2}. The CK-MB activity is obtained by multiplying the CK-B activity by two.

REAGENT COMPOSITION

R1	Imidazol, PH 6.7	125 mmol/L
	D-Glucose	25 mmol/L
R2	N-Acetyl- Cysteine	25 mmol/L
	Magnesium acetate	12.5 mmol/L
	NADP	2.52 mmol/L
	EDTA	2.02 mmol/L
	Heaokinase	≥ 6 800 U/L
	Anti – human polyclonal CK -M antibody (sheep) sufficient to inhibit up to 2000 U/L of CK-MM	
R2	ADP	15.2 mmol/L
	AMP	25 mmol/L
	di-Adenosine -5- peentaphosphate	103 mmol/L
	Glucose – 6- phosphate dehydrogenase	≥ 8 800 U/L
	Creatine phosphate	250 mmol/L

REAGENT PREPARATION

Both reagents are ready to use.

REAGENT STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm ≥ 1,2

SPECIMEN

Serum free of hemolysis or heparin plasma:

Stability 7 days at 2-8°C, protected from light

. CK-MB activity decreases a 10% after 24 hours at 4°C or 1 hour at 25°C

PRECAUTION

Working reagent (WR): Mix 4 volumes of reagent 1 with 1 volume of reagent 2. Stability: 7 days at 2-8°C or 12 hours at room temperature

(20-25°C).

PROCEDURE

1. Assay condition:

Wavelength 340 nm
Cuvette 1 cm light path
Constant temperature 25°C / 30 °C/ 37 °C

2. Adjust the instrument to zero with distilled water or air.

3. Pipette into a cuvette:

WR (ml)	1.0
Sample (μL)	40

4. Mix and incubate 10 minutes.

5. Read initial absorbance (A) of the sample, start the stopwatch and read again after 5 minutes (A2). 6. Calculate the difference between absorbances ΔA= A2 – A1.

CALCULATION

$$\Delta A \times 825 = U/L \text{ of CK - MB} \qquad \Delta A \times 1651 = \frac{U}{L} \text{ of CK - MB}$$

Calculating factor in automatic analyzers by kinetic method (ΔA/min.) is 8255.

Units: One international unit (IU) is the amount of enzyme that transforms 1 μmol of substrate per minute, in standard conditions. The concentration is expressed in units per litre of sample (U/L).

Temperature conversion factors

To correct results to other temperatures multiply by:

Assay temperature	Conversion factor to		
	25 °C	30°C	37°C
25 °C	1.00	1.53	2.38
30 °C	0.65	1.00	1.56
37 °C	0.42	0.64	1.00

NORMAL RANGE

The suspicion of myocardial damage is based on the three following factors:

	25 °C	30 °C	37 °C
CK -MB	> 10 U /L	> 15 U/L	> 24 U/L
Total CK	25 °C	30 °C	37°C
Men, up to	80 U/L	130 U/L	195 U/L
Women , up to	70 U/L	110 U/L	170 U/L

$$\frac{\text{CK - MB Activity}}{\text{CK Total Activity}} \times 100 = 6 - 25\% \text{ CK - MB Activity in the sample}$$

These values are for orientation purpose. Each laboratory should establish its own reference range.

LINEARITY

From detection limit of 1.9 U/L to linearity limit of 318 U/L .

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

QUALITY CONTROL

All control sera with CK-MB value estimated by this method can be used.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

- 1- It is important to ensure the working reagent and nitrite reagent are thoroughly mixed before adding the sample.
- 2- Bilirubin levels may be reduced if the sample is exposed to light.



Haemolytic sample will also show low value.

SYMBOL ON LABELS

Symbols	Signify	Symbols	Signify
	Catalogue Number		Pack Size
	Expiry Date		Volume
	Storage Condition		Lot Number
	Instruction for Use		In Vitro Diagnostics
	Manufacturing Date		Manufacturer
	Number of Tests		For Single Use Only
	EC Representative		European conformity

BIBLIOGRAPHY

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommendation pour la mesure de la concentration catalytique de la créatinine kinase dans la sérum humain. Ann. Biol. Clin.,40, (1482), 87.
8. Neumeier, D., Prellwitz, W., Würzburg, U. et coll. Determination of creatine kinase isoenzyme MB activity in serum using immunological inhibition of creatine kinase M subunit activity. Activity kinetics and diagnostic significance in myocardial infarction, Clin. Chim. Acta, 73, (1976), 445.

