

PRODUCT CODE
GP006

INTENDED USE

This reagent is intended for in vitro quantitative determination the amount of fructose in seminal fluid.

PRINCIPLE

Fructose reacts with boiled concentrated hydrochloric acid. This converts fructose to Oxymethyl Furfural, which then condensed with resorcinol to form a red color complex. The absorbance of this red complex is proportional to fructose concentration in the complex.

REAGENT COMPOSITION

Reagent 1

Acidic resorcinol 50 mg/dL

Reagent 2 (Standard)

Fructose 500 mg/dL

REAGENT PREPARATION

All reagents are ready to use.

REAGENT STORAGE AND STABILITY

Reagent and standard are stable up to the expiration date given on label when stored at 2-8°C

PRECAUTION

To avoid contamination, use clean laboratory wares.
Avoid direct exposure of reagent to light.

SAMPLE

Semen, centrifuge the semen sample to sediment the sperms, then take the supernatant

EQUIPMENT

Centrifuge, pipette

ASSAY

Wavelength 546nm (530 -550nm)
Cuvette 1 cm light path
Incubation Temperature Boiling Water - bath
Measurement Against reagent blank

PROCEDURE

Pipette into cuvettes	Blank	Standard	Sample
D.Water	20 µL	--	--
Reagent 2 (standard)	--	20 µL	--
Sample	--	--	20 µL
Reagent 1	1000 µL	1000 µL	1000 µL

Mix, and incubate in boiling water bath for 7-10 minutes. Cool directly all tubes under tap water, and read the absorbance of all tubes against reagent blank.

CALCULATION

$$\text{Fructose in semen. (mg/dL)} = \frac{\Delta A \text{ sample}}{\Delta A \text{ standard}} \times 500 \text{ (Std.Conc)}$$

LINEARITY

This procedure is linear up to 1200 mg/dL. For higher range, the sample should be diluted with normal saline and multiply the result by the dilution factor.

NORMAL RANGE

The normal range for fructose level of normospermic in men is 120-500 mg/dl or 6.7-27.8 mmol/L.

NOTE

Use clean disposable pipette tips for its dispensation.

SYMBOL ON LABELS

 IVD	in vitro diagnostics		manufacturing date
 LOT	lot number		expiry date
 REF	catalogue number		manufacturer
	temperature limit		instruction for use

BIBLIOGRAPHY

Text book of Medical Laboratory Technology; Praful B.Godkar
Foreman, D et al. Analytical Chem.1973.56, 584-590